

## NDV Antibody Test Kit

It is a rapid, sensitive and highly specific ELISA kit for measuring NDV antibody levels in chicken serum samples following vaccination. This ELISA test kit is used for the effective monitoring of flock health performance.

### Highlights:

- ELISA plates are precoated with inactivated antigen, blocked and stabilized. This reduces the number of steps at the time of testing.
- The entire procedure can be completed in about 2½ hours.
- The kit is provided with easy-to-use software for the interpretation of the results.

### Materials Provided:

- NDV Coated Plate
- Negative Control Serum
- Positive Control Serum
- Rabbit anti-chicken IgG HRP Conjugate
- Sample Dilution Buffer
- Stop Solution
- TMB/H<sub>2</sub>O<sub>2</sub> (20X)
- Wash Solution (10X)
- Software for the interpretation of results

### Specifications:

<b>Sample type</b>	: Chicken serum
<b>Total time</b>	: 2 hours 30 minutes
<b>Intra-assay CV percent</b>	: < 10%
<b>Inter-assay CV percent</b>	: < 10%
<b>Wavelength</b>	: 450 nm as primary filter and 650 nm as a reference filter.
<b>Stability</b>	: 6 months at 4°C

### Ordering Information:

Product	Size	Cat #
NDV Antibody Test Kit (for 450 tests)	1 Kit	105937