

GeNei™ Whole Blood DNA Extraction Teaching Kit

Description:

The kit provides the user "hands on experience" on genomic DNA extraction from fresh blood. It's a guanidium thiocyanate solution based cell lysis method, requiring no phenol-chloroform extraction.

The kit provides materials sufficient for 10 experiments


Ordering Information:

Product	Size	Cat #
GeNei™ Whole Blood DNA Extraction Teaching Kit (Consumables for 10 experiments)	1 Pack	107712
GeNei™ Whole Blood DNA Extraction Teaching Kit (Consumables for 20 experiments)	1 Pack	107713

Materials Provided:

- Solution A and B
- Solution C
- Alcohol
- Agarose
- 50X TAE
- Gel loading buffer
- 1.5 ml vials

Note: UV transilluminator, Blood and EtBr are required.



E D U C A T I O N A L P R O D U C T S