

Protein Desalting Columns


Description: These are ready to use gel filtration columns with gel of very low exclusion limit, suitable for rapid desalting or buffer exchange of proteins. These columns provide a simple and time saving alternative to dialysis for smaller sample volumes. Column of four sizes are provided to suit different volumes.

Specifications


Desalting column	105325	105326	105327	105328
Bed volume (ml)	2	5	10	20
Void volume (ml)	0.8	1.8	3.6	7.5
Sample volume (ml)	< 0.5	0.5-1	1-2.5	2-6

Removal of NaCl from BSA using Bangalore Genei's desalting columns

0.5 ml of 1% BSA in 1M NaCl
passed through DS-01 Column



5 ml of 1% BSA in 1M NaCl
passed through DS-04 Column



Ordering Information:

Catalogue Number	Product	Pack Size
Desalting Columns		
105325	Desalting Column (2 ml)	4 Nos
105326	Desalting Column (5 ml)	4 Nos
105327	Desalting Column (10 ml)	4 Nos
105328	Desalting Column (20 ml)	4 Nos
105329	Desalting Columns Pack (Contains 1 each of 2, 5, 10 & 20 ml columns)	1 Pack